

THANK YOU FOR YOUR ORDER

User Manual

Model Code: 3494

### Contents

1	Front	Cover
1	FIOH	COVE

-

2 Safety

Instructions

-

2 Warnings

-

3 Specification

\_

3 Key Functions

-

3 Operating

Instructions

-

4 Battery Safety

-

4 Back Cover

End OYPLE



# **Baby Scales**

Manual Version: **v1.0** 

Model Code: 3494

Oypla.com LLP / www.oypla.com

### Contents

### 1 Front Cover

-

2 Safety
Instructions

\_

2 Warnings

-

3 Specification

-

3 Key Functions

-

3 OperatingInstructions

\_

4 Battery Safety

-

4 Back Cover

End

# Safety Instructions

#### Notice: Ensure you read and fully understand instructions before use

While every attempt is made to ensure the highest degree of protection in all equipment, we cannot guarantee freedom from injury. The user assumes all risk of injury due to use. All merchandise is sold on this condition, which no representative of the company can waive or change.

- Only use the scales on a clean, flat table or surface.
- To ensure a long life, do not expose the scales to water or extreme heat or cold.
- This product is a high precision electronic weighing device, do not press hard or drop the scales, as this may cause damage.
- This product is not a toy, and is only to be used for its intended purpose.
- Do not carry out alterations or modifications to this product.
- Clean only using a soft, damp cloth.

- If the surface of the scales is dirty, clean the scales with a soft cloth. Do not use any chemical cleaning agents, and do not allow water to enter the body of the scales.
- If the display value is incorrect or "patchy", ensure that the batteries are adequately charged and make good contact with the terminals.
- Replace the batteries if they have been in use for a lengthy of period, and try again.
- If the product is damaged or defective, please contact <u>cservice@oypla.com</u>

# Warnings

⚠ Never leave the baby unattended whilst using the scales.

⚠ When cleaning the scales, never allow water to penetrate through into the body or battery compartment.

⚠ Be wary of the baby rolling or moving whilst on the scale – support them so that they remain still.

Manual Version: v1.0

Model Code: 3494

### Contents

1 Front Cover

2 Safety

Instructions

2 Warnings

Specification 3

3 **Key Functions** 

3 Operating Instructions

**Battery Safety** 4

Back Cover 4

End

## Specification

These scales are equipped with a high precision "strain gauge" sensor, and a highperformance CPU to ensure weighing accuracy.

Maximum Capacity: 20kg

Minimum Capacity: 10g

Resolution: 10g

Tare Function

Auto Power Off

Low Voltage Indicator

Power Source: 2 x AA Batteries

Operating Temperature Range: 10°C - 40°C

Display: LCD

Environmental Humidity: <80% R.H.

# **Key Functions**

- TARE Function If you use a towel to lay the baby on, for example, put the towel on the platform and then press the button to set the display to "0.00".
- Unit of Measurement Press this button to adjust the unit of measurement: Kilogram (kg), Pound (lb), Stone (st).
- Display This is the main display area on the scales this shows the weight of the baby.

## Operating Instructions

- Remove the scales from the packaging. Open the batter cover and install the batteries, Put the scales on a flat table, making sure it is stable.
- Press the "ON" button to turn on the power. "HI" appears on the display screen to greet you, followed by "0.00".
- Press the UNIT key to choose your desired measurement.
- Place the baby on the scales; the LCD display will show the baby's weight.
- Press the "OFF" key for 3 seconds, and switch off the scale.
  - o Note: Use two AA batteries. Open the battery cover on the back and put the two AA batteries into the battery holder, making sure the poles are aligned correctly before closing the battery cover firmly again.

Manual Version: v1.0

Model Code: 3494



### Contents

4	Front	C
	Front	( OVAr
	HOHL	

\_

2 Safety

Instructions

-

Warnings

-

3 Specification

-

3 Key Functions

-

3 Operating

Instructions

-

4 Battery Safety

-

4 Back Cover

End

# Additional Safety Instructions (Batteries)

- Use only the batteries mentioned in the manual.
- Do not install batteries in reverse polarity.
- Do not short-circuit or disassemble the batteries.
- Do not expose the batteries to water.
- Do not expose the batteries to fire or excessive heat.
- Batteries are prone to leakage when fully discharged. To avoid damage to the product, remove the batteries when leaving the product unattended for longer periods of the time
- If liquid from the batteries comes into contact with skin or clothing, immediately rinse with fresh water.

#### HAVE A QUESTION ABOUT YOUR PURCHASE?

Our dedicated customer services team are happy to help. Contact them via:

Telephone: 020 3600 22 55
Email: cservice@oypla.com
WWW.OYPLA.COM

#### CONSIDER THE ENVIRONMENT!

In the United Kingdom, approximately 5 billion corrugated boxes are used per year, amounting to around 83 per person! By recycling the cardboard box in which this product was packed, you contribute to the preservation and sustainability of the environment.

If you see a recycling logo on the packaging of your product, such as the below examples, that piece of packaging is recyclable. If there is not a recycling logo, please check with your local authorities before recycling.









PLEASE RETAIN A COPY OF THESE INSTRUCTIONS FOR FUTURE CONSULTATION



AN OFFICIAL OYPLA BRANDED PRODUCT

#### **ENVIRONMENTAL PROTECTION**

Products bearing the symbol shown below means that used electrical and electronic equipment (WEEE) should not be mixed with general household waste. For proper treatment, recovery and recycling, please take this product(s) to designated collection points where it will be accepted free of charge. Disposing of this product correctly will help save valuable resources and prevent any potential negative effects on human health and the environment, which could otherwise arise from inappropriate waste handling. Please contact your local authority for further details of your nearest designated collection point. Penalties may be applicable for incorrect disposal of this waste, in accordance with your national legislation.



Manual Version: **v1.0** 

Model Code: 3494