



Out of this world savings.

---

## User Manual

Model: 3494

---



## Infant/Pet Scales

---

CE

UK  
CA

---

# Contents

Safety Information .....	2
Operating Instructions.....	3
Technical Data.....	4

---

## Safety Information

### Please ensure you read and fully understand instructions before use

While every attempt is made to ensure the highest degree of protection in all equipment, we cannot guarantee freedom from injury. The user assumes all risk of injury due to use. All merchandise is sold on this condition, which no representative of the company can waive or change.

- 
- Only use the scales on a clean, flat table or surface.
  - To ensure a long life, do not expose the scales to water or extreme heat or cold.
  - This product is a high precision electronic weighing device, do not press hard or drop the scales, as this may cause damage.
  - This product is not a toy, and is only to be used for its intended purpose.
  - Do not carry out alterations or modifications to this product.
  - Clean only using a soft, damp cloth.
  - **Warning! Never leave the baby unattended whilst using the scales.**
  - **Warning! When cleaning the scales, never allow water to penetrate through into the body or battery compartment.**
  - If the surface of the scales is dirty, clean the scales with a soft cloth. Do not use any chemical cleaning agents, and do not allow water to enter the body of the scales.
  - If the display value is incorrect or "patchy", ensure that the batteries are adequately charged and make good contact with the terminals.
  - Replace the batteries if they have been in use for a lengthy of period, and try again.
  - If the product is damaged or defective, please contact [cservice@oypla.com](mailto:cservice@oypla.com)
  - **Warning! Be wary of the baby rolling or moving whilst on the scales – support them gently so that they remain still.**
-

---

# Operating Instructions

**The below information is designed to help you with the operation of your Oypla product**

If the product is damaged, defective or you require more information, please contact our dedicated team at [cservice@oypla.com](mailto:cservice@oypla.com).

---

## Key Functions

- TARE Function – If you use a towel to lay the baby on, for example, put the towel on the platform and then press the button to set the display to “0.00”.
- Unit of Measurement – Press this button to adjust the unit of measurement: Kilogram (kg), Pound (lb), Stone (st).
- Display – This is the main display area on the scales – this shows the weight of the baby.

## Operation

- Remove the scales from the packaging. Open the battery cover and install the batteries. Put the scales on a flat table, making sure they are stable.
- Press the “ON” button to turn on the power. “HI” appears on the display screen to greet you, followed by “0.00”.
- Press the UNIT key to choose your desired measurement.
- Place the baby on the scales; the LCD display will show the baby’s weight.
- Press the “OFF” key for 3 seconds, and switch off the scale.
  - Note: Use two AA batteries.

## Batteries

Open the battery cover on the back and put the two AA batteries into the battery holder, making sure the ends are aligned correctly before closing the battery cover firmly again.

---

## Technical Data

---

These scales are equipped with a high precision "strain gauge" sensor, and a high-performance CPU to ensure weighing accuracy.

- Maximum Capacity: 20kg
- Minimum Capacity: 10g
- Resolution: 10g
- Tare Function
- Auto Power Off
- Low Voltage Indicator
- Power Source: 2 x AA Batteries
- Operating Temperature Range: 10°C - 40°C
- Display: LCD
- Environmental Humidity: <80% R.H.

---

## Have a Question About Your Purchase?

**Our dedicated customer services team are happy to help. Contact them via:**

Telephone: 020 3600 22 55  
Email: [cservice@oypla.com](mailto:cservice@oypla.com)  
Live Chat: [www.oypla.com](http://www.oypla.com)

---

## Consider the Environment!

By recycling the cardboard box in which this product was packed, you contribute towards helping the environment. If you see a recycling logo such as the below on the packaging, it can be recycled. If a recycling logo is not present, please check with your local authority before recycling.



---

## Environmental Protection

Products bearing the symbol shown below means that used electrical and electronic equipment (WEEE) should not be mixed with general household waste. For proper treatment, recovery and recycling, please take this product(s) to designated collection points where it will be accepted free of charge. Disposing of this product correctly will help save valuable resources and prevent any potential negative effects on human health and the environment, which could otherwise arise from inappropriate waste handling. Please contact your local authority for further details of your nearest designated collection point. Penalties may be applicable for incorrect disposal of this waste, in accordance with your national legislation.



---

**PLEASE RETAIN A COPY OF THESE INSTRUCTIONS FOR FUTURE CONSULTATION**

An Official Oypla Branded Product

Imported by Oypla.com LLP, Uxbridge, UB8 2FX

**User Manual**

Model Code: 3494

[Oypla.com](http://Oypla.com) LLP / [cservice@oypla.com](mailto:cservice@oypla.com)